PHARMAC

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
STAINLESS STEEL
BIOREACTOR SYSTEM

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XBioR

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SUSPENSION SOLUTIONS
HIGHLY CUSTOMISED
AUTOMATED

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PHARMAC is committed to providing professional clean liquid system engineering services for the pharmaceutical, daily chemicals, food, synthetic biology and related industries. We have been dedicated to driving breakthroughs in fluid process technology and the development of proprietary process equipment in order to achieve high-quality and sustainable localized manufacturing of clean fluid equipment.

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We are capable of providing customers with professional fluid process equipment such as bioreactors, formulation systems, fermentation systems, PW & WFI generation & distribution systems.

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Leading the way in the clean liquid system, making contributions in the advancement of life sciences.

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Recombinant protein production in mammalian cell culture expression

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Pharmac's XBioR series stainless steel bioreactor systems, with scales ranging from laboratory-scale to industrial-scale (5L-15KL). Using CFD simulation technology, we optimize gas-liquid, liquid-liquid, and solid-liquid mixing effects, featuring low shear, high mixing intensity agitation design. We support various gas combination aeration strategies and precise temperature control, equipped with comprehensive online monitoring technologies to effectively enhance cell culture density and expression performance.

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In addition, we offer three-dimensional design of process equipment centered around bioreactors. For large-scale equipment, we can implement black-and-white partition designs to mitigate project construction risks and reduce operational costs. Our company also collaborates with renowned domestic universities to provide customers with certain cell culture media and process validation services.

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
15,000L Ultra Large Scale Stainless Steel Bioreactor

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
XBioR series has achieved:
5L-5000L standard stainless steel bioreactor system
5L-15,000L customized bioreactor system
Upstream complete equipment solution

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CQA / CPP

- | | | |
|---|--|--|
| <ul style="list-style-type: none"> ✓ pH ✓ Pressure ✓ Cell sensor (option configuration) ✓ Automatic temperature control function ✓ (Automatic) Defoaming function ✓ (Automatic / Associated) pH Control ✓ Multiple feeding strategies ✓ Prevent clogging function of exhaust filter ✓ Automatic feeding function ✓ Minimised risk of operator | <ul style="list-style-type: none"> ✓ DO ✓ Exhaust gas analyzer ✓ (Automatic) CIP/SIP function ✓ (Automatic) Pressure control function ✓ (Automatic/Associated) DO control ✓ (Automatic/Associated) Feeding Control ✓ Prevent medium reflux function ✓ Automatic seed transfer function ✓ Automatic lifting function ✓ Gas flowrate (automatic) control | <ul style="list-style-type: none"> ✓ Temperature ✓ Raman spectra(option configuration) ✓ (Automatic) CIP/SIP function ✓ (Automatic) Pressure control function ✓ (Automatic/Associated) DO control ✓ (Automatic/Associated) Feeding Control ✓ Prevent medium reflux function ✓ Automatic seed transfer function ✓ Automatic lifting function ✓ Gas flowrate (automatic) control |
|---|--|--|

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
- | | | |
|--|---|--|
| <ul style="list-style-type: none"> ✓ PLC or DCS ✓ SCADA (option) ✓ GAMP5 ✓ ISA S88 ✓ Multi level permission | <ul style="list-style-type: none"> ✓ Recipe management function ✓ Audit tracking function ✓ Automatic alarm function ✓ Real time/historical curve query function ✓ Trend curve | <ul style="list-style-type: none"> ✓ Feeding data record ✓ Batch report function |
|--|---|--|

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Black and white partition design reduces
project construction risks and operation
costs

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In terms of engineering design:

- ✓ Modular design
- ✓ Standardized design
- ✓ Redundancy design of key components
- ✓ Black and white partition design

Equipment features:

- ✓ Reasonable 3D design reduces culture media pipeline residue and production costs
- ✓ Process development, process optimization, high and low grade different configuration
- ✓ Comply with ISPE / ASME BPE and other requirements
- ✓ Dual exhaust gas filtration
- ✓ Replaceable gas sparger to satisfy customer's exploration and optimization of the production process
- ✓ Various expansion possibilities

Agitator options:

- ✓ Mixing machine
- ✓ Magnetic stirring on top
- ✓ Mechanical mixing
- ✓ Magnetic stirring

Impeller blades options:

- ✓ Elephant ear blade
- ✓ Hydrofoil blade
- ✓ Four wide blade propulsion blade, etc
- ✓ Multiple blade forms and combinations

从设计到落地, 助力客户商业化进程
From design to implementation, supporting clients in commercialization

拉曼光谱
Raman Spectra

产品视频
Product Video

我们的最优解 OUR OPTIMAL SOLUTION

标准化 STANDARDIZATION

<p>设计与制造符合cGMP指南和ASME BPE标准 The design and manufacture complies with cGMP guidelines and ASME BPE standards</p>	<p>备件的通用性 The versatility of spare parts</p>	<p>便于维护及快速安装的撬装模块化设计 Skid-mounted modular design for easy maintenance and rapid installation</p>	<p>基于ISA 88批次控制系统的自动化设计 Automation design based on ISA 88 batch control system</p>	<p>DQ/IQ/OQ/FAT/SAT符合GEP, GMP规范 DQ/IQ/OQ/FAT/SAT complies with GEP and GMP Standards</p>
--	---	---	---	---

定制化 CUSTOMIZATION

根据客户的需求定制生物反应器
Customized bioreactors according to client's requirements

<p>自动SIP(原位灭菌)/CIP(原位清洗) Automatic SIP (Sterilization in place) CIP (Cleaning in place)</p>	<p>多路进料/补料, 收获和无菌采样系统 Multi-feeding line, harvesting and aseptic sampling system</p>
<p>温度、压力、pH、DO、拉曼光谱等多种PAT技术的应用 The application of various PAT techniques, such as temperature, pressure, pH, DO, and Raman Spectra.</p>	<p>滤壳提升设计 Filter housing lifting design</p>
<p>早期设计阶段运用CFD(计算流体力学)模拟 Utilization of CFD (Computational Fluid Dynamics) simulations in the early design phase.</p>	<p>双尾气处理系统(加热/冷却) Dual exhaust gas treatment system (heating / cooling)</p>
<p>稳定的温控 Stable temperature control</p>	<p>高度灵活的进气系统 High flexible gassing strategy</p>
<p>低剪切, 强混合的搅拌系统 Low shear, high mixing intensity agitation system</p>	

XBioR 系列参数表 XBIOR SERIES PARAMETER TABLE

容器 Vessel <ul style="list-style-type: none"> 工作体积 Working volume 全容积 Total volume 最小工作容积 Minmum liquid working volume 容器高径比 Overall geometry H:D ratio 设计压力 Design pressure 设计温度 Design temperature 	<ul style="list-style-type: none"> 5L-30KL 7L-40KL ~30% 1.1:1~2:1 -1~3 barg 5-150°C
搅拌 Agitator <ul style="list-style-type: none"> 桨叶形式 Impeller type 桨叶数量 Impeller quantity 搅拌电机 Agitation motor drive 电机功率 Motor power 搅拌器类型 Agitator type 	<ul style="list-style-type: none"> 象耳桨或莱宁水翼桨 Elephant ear or lightning impellers 2 380 V AC 3/220 V AC 根据应用条件 According to the application condition 上/下机械或磁力 Top / bottom mechanical or magnetic
夹套 Jacket <ul style="list-style-type: none"> 夹套类型 Jacket type 温度控制 Temperature control 	<ul style="list-style-type: none"> 全夹套/蜂窝夹套/半管夹套 Full/ Dimple/ half-pipe jacket 冷/热板换循环控温或电加热控温 Closed loop system with recirculation pump and heat exchanger for heating and cooling or electrical heating
过程监控仪表 Process monitoring instrument <ul style="list-style-type: none"> pH Probes DO Probes MFC 泡沫电极 Foam probes 压力传感器 Pressure transmitter 液位/重量监测Level/Weight 温度监测Temp. 可选 Option 	<ul style="list-style-type: none"> 2(可备用接口) (Spare nozzle is optional) 2(可备用接口) (Spare nozzle is optional) 流量计+控制阀或MFC Flowmeter+control valves or MFC 1(可备用接口) (Spare nozzle is optional) 1(可备用接口) (Spare nozzle is optional) 差压液位计或称重模块 Differential pressure level transmitter or load cell 单支双芯(可备用接口) Single with double core (Spare nozzle is optional) CO₂ 电极/细胞密度电极/拉曼电极 CO₂ Probes / Cell density Probes / Raman Probes
通气 Gassing strategy <ul style="list-style-type: none"> 表层通气 Surface 深层通气 Deep 空气 Air 氧气 O₂ 二氧化碳 CO₂ 氮气 N₂ 过滤器 Filter 	<ul style="list-style-type: none"> 管道直通 Open pipe 环形分布器 Ring sparger 微泡管 Microsparger 按需求配置 Configure as required 按需求配置 Configure as required 按需求配置 Configure as required 按需求配置 Configure as required 0.2um 除菌过滤器 Microbial filter
尾气 Exhaust <ul style="list-style-type: none"> 换热器 Heat-exchanger 过滤器 Vent Filter 尾气分析 Vent gas analysis 	<ul style="list-style-type: none"> 排气冷却器和加热器可选配 Exhaust cooler and heater for option 0.2um 除菌过滤器(一用一备可选) Dual 0.2um microbial filter (One use, one prepare for option) 预备用接口 Spare nozzle for reservation
取样 Sampling <ul style="list-style-type: none"> 取样阀手动/自动 Sampling valve Manual/Automatic CIP/SIP 	<ul style="list-style-type: none"> 手动/自动 Manual/Automatic CIP/SIP
进/补料 Feeding <ul style="list-style-type: none"> 进料1 Feeding line 1 进料2 Feeding line 2 进料3 Feeding line 3 进料4 Feeding line 4 进料5 Feeding line 5 进料6 Feeding line 6 进料7 Feeding line 7 进料8 Feeding line 8 灌流 Perfusion 	<ul style="list-style-type: none"> 接种 Seed Feeding 培养基进料 Medium Feeding 蠕动泵/流量计可供选择 Peristaltic pumps / Flowmeter for option 蠕动泵/流量计可供选择 Peristaltic pumps / Flowmeter for option 蠕动泵/流量计/备用口可选 Peristaltic pumps/ Flowmeter / Spare nozzle for option 蠕动泵/流量计/备用口可选 Peristaltic pumps/ Flowmeter /Spare nozzle for option 可备用接口 Spare nozzle for option 可备用接口 Spare nozzle for option 可备用接口 Spare nozzle for option
控制系统 Control system <ul style="list-style-type: none"> 系统 system 显示 Display 编程 Programming 	<ul style="list-style-type: none"> Simens Wincc or PCS7+batch 触摸屏或电脑 Touch screen or IPC ISA-88 GAMP5

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- 制造始于1995年
Our history of manufacturing can be traced back to 1995
- 20年+, 洁净容器/换热器/模块设计及加工制造历史
20+ years' experience in clean vessel/heat exchanger/module design and manufacturing
- 容器及模块制造工艺流程的过程控制
Process control of vessel and module manufacturing process flow
- 国内少有可生产钛材/哈氏合金等特材压力容器的厂家
One of the few domestic manufacturers with the capacity to produce pressure vessels made of special materials such as titanium and Hastelloy

认证 AUTHENTICATION

- 压力容器特种设备生产许可证
Pressure Vessel as Special Equipment Production License
- ASME (美国机械工程师协会)U钢印授权证书
"U" Steel Seal Certificate of Authorization from ASME
- 压力设备(PED) 认证
Pressure Equipment Directive (PED) Certification
- 可提供满足 GMP、CFDA、FDA、EMA、WHO等体系的文件
We can provide documents that meet the requirements of GMP, CFDA, FDA, EMA, WHO and other systems

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Pharma Champ Pharmaceutical Engineering (ShangHai) Co., LTD.

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10F, Bldg.3A Tianhua Information Park, No.299, Longcao Road, Xuhui District, Shanghai

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Haixin Sanitary Technology Engineering (Guangzhou) Co., LTD.

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Production Site : No.7 Meide Er Rd, Zhujiang Industrial Zone, Nansha District, Guangzhou

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Production Site : No. 46, Shijing Road, Haichang Street, Haining City, Jiaxing City, Zhejiang Province

PHARMAC EMEA

William Craig

Polizeigasse 4, 86720 Noerdlingen, Bayern, Germany

+49 (0)160 894 8284

info@pharmac-emea.com

PHARMAC ASIA

21 Bukit Batok Crescent #09-79, Wcega Tower, Singapore

www.pharmachamp.com

PHARMAC-2420-CE800-20241117

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Pharma Champ Pharmaceutical Equipment (Zhejiang) Co., LTD.